# ZQuiet® S.A.

## Oral Appliance System for Treating Snoring & Obstructive Sleep Apnea

### Patient Instructions for Use and Care

#### Device Description

The ZQuiet® S.A. is a pre-fabricated, non-custom oral appliance that provides advancement of the lower jaw to treat snoring and Obstructive Sleep Apnea (OSA). The ZQuiet® S.A. comes in four sizes that provide progressive advancement of the jaw in 2mm increments. Each mouthpiece consists of an upper and lower tray that fits comfortably over the teeth and is joined by a flexible hinge to allow the mouth to open and close. The trays are designed to position the lower jaw slightly forward to open the airway in the back of the throat.

#### Prescription Device

**Caution:** The sale of this device is by prescription only.

#### Indications for Use

ZQuiet® S.A. is intended for the treatment of nighttime snoring and mild to moderate obstructive sleep apnea in adults.

#### Instructions

1. Follow your physician's instruction for using the ZQuiet® S.A. device.
2. Prior to first use, thoroughly clean ZQuiet® S.A. according to the cleaning instructions in the following section.
3. There are 4 mouthpiece sizes (marked A, B, C & D), each providing an increasing amount of jaw advancement. Device A has the least advancement and Device D has the greatest. Each are identified with both a number and a corresponding letter (0&A, 2&B, 4&C, and 6&D) inside of a circle on the under side of the upper tray (Diagram 1). Unless otherwise instructed, select the mouthpiece with the least advancement marked with both the "0" and "A" as the initial mouthpiece to begin treatment.

#### What to expect from your new ZQuiet® S.A.

Just like anything new, for many people ZQuiet® S.A. takes a little time to get used to wearing. After all, most of us are not accustomed to sleeping with a mouthpiece. So it may take some users up to a few nights to get over any possible drooling, unusual sensations in the teeth or jaw, or even slight discomfort. ZQuiet® S.A. is designed to be sleek and low profile so it should fit comfortably in your mouth during sleep. Stick with it and you will be rewarded with quiet, restful sleep. For any jaw discomfort you may experience when you wake up, or other concerns, please review the "Tips" section.

#### Possible Side Effects

Some of the following possible side effects can be associated with use of The ZQuiet® S.A. Many of these are normal during the adjustment period with any oral appliance and should disappear with continued use:

- Slight tooth or gingival discomfort due to pressure from the appliance.
- Increased salivation initially.
- Inadvertent removal of the appliance during sleep.
- Slight jaw soreness or tightness.
- A sensation of bite change after removing the device which is normal and subsides within an hour or two.

#### Cleaning Instructions

Clean your ZQuiet® S.A. each morning to minimize harmful germs and bacteria as follows:

1. Fill a cup, bowl or basin with very warm water and a small amount of mild soap or detergent.

4. Before every use, rinse thoroughly with warm water. When inserting ZQuiet® S.A., hold the mouthpiece in the folded position with the upper tray with the "Contoured Cut-Out" facing up. Insert ZQuiet® S.A. with the hinges towards the rear of your mouth (see Diagram 2) until your teeth settle into the upper and lower trays. You should notice that your lower jaw will feel slightly more forward than normal.

5. If snoring persists, it is indicative that more jaw advancement may be necessary. Proceed to use the device that provides the next level of advancement. If snoring still persists, repeat this process again for the next level of advancement until satisfactory results are achieved. If you experience jaw discomfort with the increased advancements, return to the previous device of lesser jaw advancement. If more serious pain persists (lasting more than 2 hours), discontinue use immediately and contact your prescribing physician or dentist.

6. Upon waking, gently remove the device and follow the recommended cleaning and storage instructions.

All patients who have been diagnosed with OSA should undergo some type of objective sleep testing to evaluate the efficacy of the treatment.

#### Contraindications

ZQuiet® S.A. is contraindicated for patients with loose teeth, orthodontia, full dentures or advanced periodontal disease. This device is also contraindicated for patients that have central sleep apnea, severe respiratory disorders, or have been diagnosed with temporomandibular joint disorders (TMD). ZQuiet® S.A. is not for use by anyone under the age of 18.

2. Using a toothbrush, scrub all surfaces of the device vigorously for about 1-2 minutes.
3. Rinse thoroughly with clean warm water to remove all soap.
4. After rinsing, shake off all excess water and pat dry with a soft clean towel.
5. Store device in the supplied case at room temperature to prevent damage.

**Important:** Do not soak ZQuiet® S.A. in mouthwash, especially if the mouthwash contains alcohol.

#### Comfort Tips

- Get used to your ZQuiet® S.A. initially by wearing it for short periods of time (30 minutes while watching TV or reading).
- Upon waking in the morning, your bite may feel slightly out of alignment. This sensation should go away during the next couple of hours; gently chewing some sugarless gum can help speed up the process.
- If you experience morning jaw soreness, use a gentle fingertip massage by placing your index and middle fingers on both sides of your jaw just in front of your ears. Bite down to feel this muscle contract so you can pinpoint exactly where to massage. Now, relax the jaw muscle and use your fingers to gently massage this muscle in a circular motion.
- If your jaw muscles are particularly sore, take a night off. Give your jaw muscles a chance to rest as you slowly become accustomed to using your ZQuiet® S.A.
- Place a warm washcloth over the affected jaw muscles for up to 20 minutes, 3 or 4 times per day.
- To provide further relief from any jaw soreness, consider

Diagram 1



Diagram 2  
Proper Insertion

#### Precautions

Each patient's medical history including history of asthma, breathing, or respiratory disorders, or other relevant health problems should be considered before prescribing the device. Proper screening and diagnosis by a qualified medical professional should be completed before prescribing the device for obstructive sleep apnea.

#### Warning

Usage of ZQuiet® S.A. may cause:

- Tooth movement or changes in dental occlusion.
- Gingival or dental soreness.
- Pain or soreness to the temporomandibular joint.
- Obstruction of oral breathing.
- Excessive salivation.

taking an over-the counter anti-inflammatory pain medication. It is important to use over-the-counter medicines responsibly and in accordance with labeling, and speak with a doctor or pharmacist if you have any questions.

- For additional comfort, you can trim your ZQuiet® S.A. with small scissors to improve its fit. After trimming and to avoid irritation to gums, use fingernail file or emery board to smooth any rough edges.
- Discontinue use and consult your doctor or dentist for acute or persistent pain lasting more than 2 hours, clicking or popping sensation in the jaws, or sudden limited jaw mobility.

#### Distributed by:

Kinwood Healthcare Limited  
10/F, Sino Cheer Plaza, 23 Jordan Road,  
Kowloon, Hong Kong  
<https://www.sleepkinwood.com/>

Copyright ©2016 Sleeping Well, LLC. All Rights Reserved  
ZQuiet® is a Registered Trademark of Sleeping Well, LLC.  
U.S. Patent #8166976, Foreign Patents Pending  
ZQSA.d-IFU 01.16

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